

PRINT DATE: Monday, 12 October 2015



► VIEW DATASHEET ON-LINE www.chinesport.com/go/01933/

01933 SHOWER MATTRESS


It allows bed-ridden, elderly or disabled patients to be washed and cared for in their personal hygiene with running water directly on their bed, without any need for transfer/lift.

Inside the "Shower Mattress", the patient can be showered with water coming from the tap via an extension hose, or alternatively, using the special pressurized water tank provided with flexible hose and shower head. Once finished, the water is drained into a suitable container.

Complete with plug and flexible hose for quick water drain.



01933 - SHOWER MATTRESS



Technical specifications and dimensions

Dimensions : 80 x 180 x 20h cm

Documentation available online

DOCUMENTS IN ITALIANO

User and Maintenance Manual - Istruzioni per l'uso  .pdf 3.2 MB

DOCUMENTS IN SPANISH

User and Maintenance Manual - Manual de instrucciones  .pdf 3.03 MB

DOCUMENTS IN ENGLISH

User and Maintenance Manual - User manual  .pdf 5.71 MB

Brochure - mod.14-2013.07.01-16.26-001-gb.pdf .pdf 622.76 KB

 The document is available for download in a reserved area

Accessories



01939 CARRY BAG

A container for holding and carrying the mattress



01940 WATER TANK

Pressurized water tank complete with flexible hose and shower head operated by button.



XBC001 SHOWER PANEL


Shower panel consisting of a stainless steel support and ABS outer panel. It has vacuum breaker valves to empty the water system completely after every use, avoiding any stagnation of water or sanitizing solutions.




01938 SHOWER EXTENSION HOSE

Flexible, 6m hose complete with button-operated shower head and tap adapter

More pictures of product



Turn the patient on his back again and place the towel in the middle of the bed



The washing operation